

# Circulating Rare Cells (CRCs) Technology



INSTRUCTION FOR USE OF ScreenCell® Molecular Biology (cat# MB 4LC)

## DNA and RNA extraction from whole blood CRCs ScreenCell® Live Cell (SC LC) Dilution Buffer (DB) preparation

Add B2 and C reagents to 150 ml of buffer D. Invert the bottle twice. The reconstituted buffer is sterilized by autoclaving and filtrated through a 0.22 µm Millipore filter unit. The reconstituted buffer must be stored at 4°C and used within 12 hours.

### Dilution of Blood with the Dilution buffer

- Collect the blood sample in an EDTA tube. Invert the tube 3 times and keep at + 4°C. Blood must be filtered within 3- 4 hours after drawing.
- Transfer 1 ml of blood into one 15 ml RNase-/DNase- free conical tube.
- Add 7 ml of LC DB.
- Close the tube, homogenize once by inverting the tube and incubate 2 minutes at room temperature.

### Blood Filtration

- Remove the protective membrane at the extremity of **module B** (1).
- Put the SC MB-LC Unit on a dry and clean place.
- Fill **module A** with the prepared 8 ml of diluted blood (2).
- Insert **module C** and push fully upward into **module B** until filtration starts (3).
- Before filtration is completed, add 1ml RNase-/DNase- free PBS 1X in module A and keep filtering to rinse out the filter until all liquid is gone.

**Note:** filtration is usually completed within approximately 50 seconds (a sample must be considered coagulated when filtration exceeds 60 seconds).

- At the end of filtration, unscrew, remove and discard **modules B including C** from **module A** (4).
- Hold the **module A** onto a Nuclease-free Eppendorf tube (5), release and insert the capsule-filter (CF) into the Eppendorf tube by pushing down a rod located at the bottom part of **module A** (6).
- Once the CF has been inserted into the tube, remove and discard the rest of **module A** (7).

### Cell lysis for DNA/RNA extraction


- In each CF, add a volume\*\* of lysis Buffer\*\*\*. Close the CF with the Eppendorf tube cap (8).
- Incubate at appropriate temperature\*\*\*.
- Centrifuge the Eppendorf tube/CF 1 min at 12 000 g (10 500 rpm) (9).
- Discard the CF, save and store the flow-through in the closed Eppendorf tube with the cap lock.

\*1.5 ml Eppendorf tube PCR Clean, cat# 3810X.

\*\*210 µL for DNA or 165 µL for RNA

\*\*\*QIAGEN: cat# 74004; 74404; 56304; 56504.

•Ordering/Technical Information: •Phone + 33 (0) 1 53 61 91 33 •Fax + 33 (0) 1 53 61 91 36



DNA/RNA PROTOCOL FOR LIVE CRCs